


Quick Start Guide

Say hello to the revolutionary
power of rechargeable
hearing aid batteries.




ZPOWER®

QUICK START GUIDE—Instructions




Plug the power cord into the USB wall adapter and into a power outlet. The charger should flash green once to indicate power is flowing to the charger.



Place your hearing aids in the designated slots in the charger—the hearing aid wires should be facing the same side of the charger as the power cable.


QUICK START GUIDE—Instructions

3



Green blinking lights will indicate the hearing aids are charging. When the lights stop blinking and turn solid green, the hearing aids are completely charged. Charging can take up to 7 hours.

4



After the lights turn solid green, your hearing aids are ready to use. Hearing aids will power on automatically and are ready to wear.

QUICK TIPS

- 1 Fully charge your hearing aids every night.
- 2 If you stop wearing your hearing aids during the day, place them back in the charger. Make sure the charger is plugged in.
- 3 Regularly use a soft tissue or cloth to keep contacts (in charger wells and on your hearing aid doors) dry and clean.
- 4 If you hear the hearing aid low battery indicator, place hearing aids in designated slots on the charger or temporarily change out the ZPower battery and replace with a zinc-air disposable. Never try to squeeze more battery life by opening and closing the battery door.
- 5 If your hearing aids are not being worn for two weeks at a time, remove the batteries. Batteries should be stored in a safe and dry location where they will not touch each other or other metal objects.
- 6 If one or both of the charging lights are solid red, it means that the charger does not recognize the battery inside the hearing aids; verify the battery is a ZPower silver-zinc battery. If it is and you're still having trouble, try removing and reinserting the battery. If the light continues to be red, consult your hearing care provider as the battery may need to be replaced. Batteries should be replaced at least once a year.



Visit www.zpowerbattery.com for a copy of the user guide, FAQs, and more helpful tips.

9-07-00112-1 A

ZPOWER®